




Fold it up and take the Mini Focus into any situation requiring efficient, high-quality ultrasound.



Mini Focus Features

- Portable color ultrasound
- High image quality
- Up to 12 MHz frequency
- 2 transducer connectors
- Transducers for: transrectal and transvaginal scanning; breast, orthopedic and musculoskeletal scanning; vascular and small part scanning; pediatric and intraoperative scanning
- Easy-to-use puncture guides that are compatible with modern sterilization techniques
- Efficient internal archiving system to save images; USB port to save images to a portable hard drive and send them via email; email images directly from the scanner*
- Options for:
 - 3D ultrasound
 - DICOM® Store and DICOM® Print
 - OB/Gyn extended measurement package



Transducer Overview


	8665	8660	8803	8667	8806	8670	8815	8817
Abdominal								
Breast								
Intraoperative								
Musculoskeletal								
Obstetrics								
Pediatrics								
Peripheral vascular								
Small part								
Transrectal								
Transvaginal								
Cardiac								
Harmonic imaging								
Puncture								
MFI								
Frequency range	2.5-5	5-8	3-5	5-10	4.3-9	5-12	5-10	2-4.3



The Mini Focus trolley (optional accessory) is compact and easy-to-manuever.

*subject to local regulation; requires SMTP mail server

B-K Medical Worldwide



Sales Offices

B-K Medical Systems, Inc.
250 Andover Street
Wilmington, MA 01887, USA
Tel: (1) 978-988-1078 / 800 876-7226
Fax: (1) 978-988-1478

B-K Medical Benelux NV
Generaal de Wittelaan 19/7
2800 Mechelen, Belgium
Tel: (32) 1528 0600
Fax: (32) 1529 1025

B-K Medical Medizinische Systeme GmbH
Pascalkehe 13
25451 Quickborn, Germany
Tel: (49) 4106 99 55-0
Fax: (49) 4106 99 55-99

B-K Medicale S.r.l.
Via Tulipani, 3/B
20090 Pieve Emanuele MI, Italy
Tel: (39) 02 90 78 13 47
Fax: (39) 02 90 78 19 05

B-K Medical Norge
Nesbru Senter, Fekjan 7B
1394 Nesbru, Norway
Tel: (47) 66 77 45 77
Fax: (47) 66 77 45 78

B-K Medical AB
Mänskärsvägen 9, 141 75
KUNGENS KURVA, Sweden
Tel: (46) 8-744 02 11
Fax: (46) 8-744 02 12

B-K Medical (UK)
11 Grove Park
Waltham Road, White Waltham
Berkshire SL6 3LW, United Kingdom
Tel: (44) 1628 825-770
Fax: (44) 1628 826-970

Your local representative

Represented in

Argentina • Australia • Austria • Bangladesh • Belgium • Bolivia • Brazil • Brunei • Bulgaria • Cambodia • Canada • Chile • P. R. of China • Colombia • Costa Rica • Croatia • Czech Republic
Denmark • Ecuador • Egypt • Estonia • Finland • France • Germany • Ghana • Greece • Guatemala • Hong Kong • Hungary • Iceland • India • Indonesia • Iran • Iraq • Ireland • Israel • Italy • Japan
Jordan • Kazakhstan • Korea • Kuwait • Laos • Latvia • Lebanon • Lithuania • Luxembourg • Malaysia • Malta • Mexico • Myanmar • The Netherlands • New Zealand • Norway • Pakistan • Peru
Philippines • Poland • Portugal • Romania • Russia • Saudi Arabia • Serbia & Montenegro • Singapore • Slovak Republic • Slovenia • South Africa • Spain • Sweden • Switzerland • Syria • Taiwan
Thailand • Tunisia • Turkey • Ukraine • UAE • UK • Uruguay • USA • Venezuela • Vietnam

World Headquarters
Mileparken 34
DK-2730 Herlev
Denmark

Tel: (45) 44 52 81 00 • Fax: (45) 44 52 81 99 • www.bkmed.com



Mini Focus 1402

Gain more opportunities with portable color




ultrasound
we specialize because you do



A new level of versatility with portable color

B-K Medical transducers come with convenient puncture guides that are compatible with modern sterilization techniques.



Mini Focus is packed with all the tools and transducers you need for both specific and general ultrasound applications. It's a compact, streamlined scanner you can fold up and take into any situation requiring fast, user-friendly ultrasound.

- Includes convex, linear and endo transducers.
- Works with a high-frequency linear transducer, which is ideal for breast and orthopedic scanning, rheumatology, sports medicine and other applications where high resolution is a demand.
- Works with a compact and easy-to-hold cardiac and emergency transducer.
- Gives you the crisp image details, deep views, and convenient puncture guides you need to gain even more confidence for general ultrasound applications.

In both general and emergency medicine, the ease-of-use and compact size of the Mini Focus are key advantages; carry it, use it quickly and then move it away when you're done.

- Easily and quickly get the clear, high-quality images you need with few adjustments, so you can save time and see more patients.
- The color Doppler feature of Mini Focus, combined with our excellent abdominal probe, lets you clearly see pooling of blood and blunt trauma to kidneys, liver and spleen.
- No need to bring in extra monitors to an emergency situation – you see everything you need to see directly on the Mini Focus, with its excellent, flat panel monitor.




Save and share information easily



The Mini Focus offers a new level of convenience in image storage and retrieval. With its internal archiving system, you can easily save ultrasound images no matter where you are and in what situations you use it.

- Save images to an internal hard drive, at the touch of a button.
- Keep track of images from the same patient over several visits, and if need be, compare later scans with earlier ones.
- DICOM® Store and DICOM® Print.
- Save images to a portable hard drive via a USB port.
- Share information with colleagues: email pictures directly from the Mini Focus*.

*subject to local regulation; requires SMTP mail server



Expand your practice capabilities

With the Mini Focus, you can expand into new application areas for ultrasound at a cost-effective price.


- A portable color system that offers new diagnostic possibilities (such as color Doppler) to perform breast, musculoskeletal, kidney, testis, penile Doppler scans, and many more.
- Includes a linear probe which is excellent for breast scanning.
- Includes the B-K Medical dedicated endfire prostate probe; developed with urologists, to give greater confidence in prostate examinations.
- Convenient biopsy system designed to make biopsies easier and quicker.



Bottom Left: Cardiac scan. Apical long axis, mitral inflow.

Bottom Center: Increased vascularity around prostate cancer.

Bottom Right: Breast scan. Normal breast.



Top: Transverse scan of the wrist.

Center: Abdominal scan of the liver using tissue harmonic imaging.

Bottom: Color Doppler showing normal flow in the kidney.

