


MASTERPULS® MP50 – The **COMPACT**
radial shock wave therapy system



MASTERPULS® MP50 – Radial shock wave therapy system with D-ACTOR® technology



Myofascial trigger points (MTrP) are the most frequent causes of acute and chronic pain in the locomotor system. Local improvement of perfusion, reduction of nociceptor activity and neuromuscular tone and local inter- and intramuscular stretching are the physiological and mechanical axioms of MTrP therapy with the MASTERPULS® MP50.



The highlights of this exceptionally efficient shock wave system include: compact design, integral compressed air supply, low capital expenditure and maintenance costs and maximum mobility thanks to the built-in compressor.

Facts and figures:

- ▶ Compact design: 9 kg, 34 x 34 x 15.6 cm
- ▶ Built-in compressor
- ▶ Shock frequency: 1 – 11 Hz
- ▶ Application pressure: 1 – 3 bar/11 MPa
- ▶ Control via LED display
- ▶ Therapy recommendations: 3 – 6 sessions at weekly intervals



- ▶ Therapy of myofascial pain
- ▶ Activation and loosening of muscle and connective tissue
- ▶ Disorders of tendon insertions and ligaments
- ▶ Acupuncture shock wave therapy